



Smart Technology, Made to Care

## Pill

Pivotell Advance

Registering  
Guide v1.0



The Pivotell Advance Pill Dispenser allows medication to be dispensed safely and securely at the exact time it is required. When the medication is ready to be taken, the device will rotate to release the medication, and will produce an audible and visual alert.

The pills can then simply be tipped from the device, to be consumed by the user. For added security, the dispenser is lockable and has a self-locking shutter, so the medication can only be accessed when it is needed, avoiding the risk of double-dosage.

The pill dispenser can be used on its own, or connected to a Notifier or Lifeline home unit, to alert if the medication has not been taken.

### Registering the Pill Dispenser to a Notifier/Lifeline

1. Open the lid of the battery compartment on the underside of the dispenser.
2. Ensure the dispenser is activated by removing the battery isolating strip.
3. Enter registration mode on the Notifier/Lifeline
4. A switch will also be found within the battery compartment.
5. Press the switch, the dispenser will register with the Notifier/Lifeline.
6. Press the switch again to ensure it activates the Notifier/Lifeline.
7. Replace the battery compartment lid, ensuring that the switch is lying across the middle of the four batteries.

Please refer to the enclosed instructions on how to configure the Pill Dispenser.