



Smart Technology,
Made to Care



CairWear Quick Reference Guide

we-cair.com

Contents

- Page 3 - How it Works
- Page 4 - Wearing Options
- Page 5 - Switching On/Transport Mode
- Page 6 - Using the Device
- Page 7 - Disabling/Enabling the Haptic Feedback
- Page 8 - Registering to a Notifier/Other System
- Page 9 - Changing the Battery



CairWear

How it Works

The CairWear is a sleek, modern pendant designed to be worn comfortably on the wrist or around the neck. It provides users with peace of mind by allowing them to easily call for help with the press of a button.

The device features haptic feedback, ensuring users receive clear confirmation that their call for help has been successfully activated. For those with dexterity challenges, an optional Easy-Press Adapter is available to make operation even simpler.

Built to withstand daily wear, the CairWear is waterproof to IP67 standards, includes a replaceable battery, and comes equipped with built-in Auto Low Battery functionality.



Smart Technology,
Made to Care

CairWear

Wearing Options

The CairWear can be worn in two ways - on the wrist or around the neck.

If the neck cord is to be used, thread the looped end of the cord through one of the designated openings on the CairWear device. Pull the rest of the cord through the loop to secure it. Once attached, the device is ready to be worn around the neck.

If the wrist strap is preferred, first detach the neck cord. Then, thread the strap through both openings on the CairWear device. Once secured, the device is ready to be worn on the user's wrist.

If the Easy-Press Adapter Jacket is required, stretch the jacket and fit it securely over the CairWear device. Once in place, the device can be worn around the neck or on the wrist, depending on the user's preference.



Smart Technology,
Made to Care

CairWear

Switching On/Transport Mode

Switching the CairWear On

To switch on the device, press the small button once. The LED will flash red, green then orange to indicate that the unit is now on and ready for use.

Entering Transport Mode

To place the CairWear into Transport Mode, press and hold the small button until the LED turns off (green for a good battery, red if the battery is low). Without releasing the small button, press and hold the manual button until the orange LED stops flashing. Release both buttons as soon as the orange LED stops flashing. The device is now in Transport Mode.



Smart Technology,
Made to Care

CairWear

Using the Device

To activate a call for help, ensure the CairWear device is switched on. Press the large button once to send a radio signal. The device will provide vibratory and LED feedback to confirm that the call has been raised. A green LED will indicate the battery is good and a red LED will indicate the battery is low.

Note: The haptic feedback (vibration) can be disabled

A battery check can also be performed by pressing the small button when the device is switched on. Pressing it once will cause the LED to display green if the battery is good or red if the battery is poor.



Smart Technology,
Made to Care

CairWear

Disabling/Enabling the Haptic Feedback (vibration)

To toggle the haptic feedback (vibration) setting between on or off, ensure the device is in Transport Mode, then press the small button twice in quick succession. The device will flash red, green, orange, then red to indicate that the setting has changed.

To confirm whether haptic feedback is enabled or disabled, place the device into Transport Mode, then switch it back on. Press the manual button. If haptic feedback is enabled, the device will flash and vibrate. If haptic feedback is disabled, the device will flash but will not vibrate.



Smart Technology,
Made to Care

CairWear

Registering to a Notifier/Other System

- Ensure the CairWear is switched on.
- Prepare the Notifier/other system by entering registration mode.
- Activate the CairWear by pressing the large button.
- The Notifier/other system should acknowledge the new device.



Smart Technology,
Made to Care

Thank you



Smart Technology,
Made to Care

www.we-cair.com

01422 399155 | hello@we-cair.com